

Product ID: 15 | Category: PC

# Laboratory Chemical Container

High-impact PC bottle, chemical-resistant, for demanding lab environments.

## General Information

MATERIAL	WEIGHT (1L)
PC (Polycarbonate)	180-250g
TEMPERATURE RANGE	NECK FINISH
-100°C to 135°C	GL45 Thread
TRANSPARENCY	RECYCLING CODE
89-91%	7 (PC)
AUTOCLAVABLE	CHEMICAL RESISTANCE
Yes (121°C, 20 min)	Excellent

## Available Specifications

Volume	Available Colors	Price Range (USD)
500ml	Clear, Blue, Green	\$4.00 - \$4.05
1L	Clear, Blue, Green	\$5.00 - \$5.05
1.5L	Clear, Blue, Green	\$6.50 - \$6.55
2L	Clear, Blue, Green	\$9.00 - \$9.05

## Key Features

- ✓ Excellent chemical resistance to acids, bases, and solvents
- ✓ High impact strength (250x stronger than glass)
- ✓ Autoclavable for sterilization (121°C, 20 minutes)
- ✓ Wide temperature range for laboratory use
- ✓ Graduation markings for accurate measurement
- ✓ Leak-proof screw cap with PTFE seal
- ✓ Crystal clear for content visibility

## Applications

Ideal for laboratory chemicals, reagents, acids, bases, solvents, buffers, and sterile solutions. Suitable for pharmaceutical, biotech, research, and industrial laboratories.

## Technical Specifications

Specification	Value
Density	1.20 g/cm³
Impact Resistance	250x glass strength
Heat Resistance	Up to 135°C
Chemical Resistance	Excellent
Wall Thickness	2.0-3.0mm
Sterilization	Steam autoclave compatible

## Chemical Compatibility Chart

Chemical Type	Compatibility	Notes
Acids (dilute)	Excellent	Up to 50% concentration
Bases (dilute)	Excellent	Up to 40% concentration
Alcohols	Good	Most common alcohols
Solvents	Good to Excellent	Varies by solvent type
Saline Solutions	Excellent	All concentrations
Buffers	Excellent	pH 2-12 range

## Quality Standards & Certifications

ISO 9001:2015  
FDA Approved  
USP Class VI  
REACH Compliant  
RoHS Compliant  
UN Approved

**Warning:** Always verify chemical compatibility for specific applications. Contact technical support for detailed compatibility data.

## Sterilization & Cleaning

- ✓ Autoclavable: 121°C, 20 minutes, 15 psi
- ✓ Dry heat sterilization: Up to 135°C
- ✓ Chemical sterilization compatible
- ✓ Wash with laboratory-grade detergents
- ✓ Rinse with deionized water
- ✓ Dry in clean environment

## Storage & Handling

- ✓ Store in cool, dry, well-ventilated area
- ✓ Keep away from incompatible materials
- ✓ Use appropriate PPE when handling
- ✓ Follow laboratory safety protocols
- ✓ Label containers clearly
- ✓ Store in chemical-resistant cabinets

## Environmental Information

### RECYCLABILITY

Yes (Check local programs)

### REUSABILITY

Multiple sterilization cycles

### RECYCLING CODE

7 (PC)

### LIFE CYCLE

3-5 years laboratory use